

KYZ2002-5L 5 Liter Oxygen Concentrator

PARAMETERS

Model Number	KYZ2002-5L
Type	Medical Oxygen Concentrator
Max.Flow	5L/min
Oxygen Density	93 ₋ 3%(V/V)
Output Pressure	20kpa ~ 60kpa
Power Consumptions	≤ 480W
Noise	≤52dB (A)
Continuous Working Time	15,000 Hours
Product Size	39*39*69cm,1 pcs/ctn
Gross Weight	18Kgs
Display	LCD